



# coDiagnostiX® Version 9.15

## LICENSE FEATURES AT A GLANCE

	PRODUCER	CLIENT	EASY	EASY CHAIRSIDE
<b>Planning</b>				
EASY user interface	●	●	●	●
Panoramic curve	●	●	●	●
Patient coordinate system	●	●	●	●
Planning of implants, abutments & sleeves	●	●	●	●
Comprehensive implant library	●	●	●	●
<b>Safety</b>				
Nerve detection	●	●	●	●
Distance warning	●	●	●	●
<b>Data Import</b>				
DICOM import	●	●	●	●
Model & intraoral scan import	●	●	●	●
CARES® / DWOS® Connect	●	●	●	●
3Shape Communicate	●	●	●	●
<b>Data Output</b>				
Printouts	●	●	●	●
iChiropro	●	●	●	●
iPad viewer	●	●	●	●
External navigation system (Minivident)	●			
<b>Surgical Guides</b>				
Surgical guide design	●	●	●	●
Surgical guide export	●			●
<b>Collaboration</b>				
Transfer cases via caseXchange™	●	●	●	
Transfer cases via archive exchange	●	●		
Multi-seat planning (network database)*	●	●		
<b>Prosthetics-driven Planning</b>				
Virtual teeth	●	●		
CARES® / DWOS® Synergy	●	●		
Third-party CAD interface	●	●		
<b>Advanced Planning</b>				
Distance measurements	●	●	●	●
Continuous distance measurements	●	●		
Angle measurements	●	●		
EXPERT user interface	●	●		
Annotations	●	●		
Augmentation	●	●		
Image management	●	●		
Implant designer	●	●		
Custom sleeve system	●	●		
Dual scan procedure	●	●		
Support for edentulous cases	●	●		
Segmentation	●	●		
Treatment evaluation*	●	●		
Virtual planning export	●			