

## PC specifications to run DWOS software

Version 2.3 4/4/2018

---

### Recommended PC specifications to run DWOS software:

- CPU: Quad Core 3.30 GHz [Intel Core i5 3.30 GHz]
- Memory RAM: 16 Gigabytes
- Graphics card: ATI Radeon HD6570, 1GB
- Hard drive (HD): 500 gigabytes 7200 rpm
- Network with Internet connection
- Monitor 22" wide 1080p (Recommended)
- USB keyboard
- USB mouse (two buttons and wheel)

### Optional components:

- 3D-mouse (3DConnexion brand - Space Navigator, Spacemouse Pro and Space Explorer compatible)

### Operating system:

- Windows 7 (64-bit edition)

### Additional specifications if the PC is driving a 5Series scanner:

- Four (4) USB ports. Two (2) of them must be free once keyboard and mouse are connected.
- One (1) 1394A firewire connector with Texas Instrument chipset (Vantek UGT-FW210 guaranteed compatible) **if and only if the 5Series is PGR firewire.**