

PC specifications to run DWOS software

Version 2.4 3/24/2020

Recommended PC specifications to run DWOS® Software

- CPU: Quad Core 3.30 GHz [Intel Core i5 3.30 GHz]
- Memory RAM: 16 Gigabytes
- Graphic card: ATI Radeon HD6570, 1GB
- Hard drive (HD): 500 gigabytes 7200 rpm
- Network with internet connection
- Monitor 22" wide 1080p (recommended)
- USB keyboard
- USB mouse (two buttons and a wheel)

Operating system

- Windows 7 (64-bit edition)
- Windows 10 (64-bit edition)

Optional components

- 3Dconnexion navigation devices

Additional specifications if the PC is driving a 5Series

- Four (4) USB ports; two (2) of them must be free once keyboard and mouse are connected.
- One (1) 1394A firewire connector with Texas Instrument chipset (Vantek UGT-FW210 guaranteed compatible) **if and only if the 5Series is PGR firewire**