

7SERIES

Model & Impression Scanner


Featuring high scanning accuracy, significant ergonomic and functional advances for important improvements in productivity and efficiency. Delivered with several DWOS design applications suited for medium to large sized dental labs. Plus, you can add more modules as your lab and requirements evolve.



Accessories included:

- ▶ An integrated on-board computer with 64-bit processor
- ▶ Impression holder
- ▶ Calibration kit
- ▶ Automatic multi-die holder (up to 30 elements scanned and designed in 12 minutes, unattended)

Choose between two DWOS packs: **Synergy** and **Productivity**



Technical specifications	7Series
Scanning volume	140 x 140 x 140 mm [5.51 x 5.51 x 5.51"]
Optical technology	Class 1 laser product, 2 high speed measuring cameras, 1 color video camera
Number of axes	5 (3 rotative 2 linear)
Embedded computer & OS	Core i7, 32 GB memory, 4 GB of dedicated RAM graphic card, Windows 10, 64 bits
Accuracy	15 microns ¹
Output format	STL & supported partner formats
Scannable materials	Plaster, wax-up and impression

¹ Based on Dental Wings testing standards



International Headquarters

Dental Wings Inc.
160 Rue St-Viateur E
Montréal QC H2T 1A8 CANADA

DENTAL WINGS and/or other
trademarks and registered trademarks
from Dental Wings mentioned herein
are the trademarks or registered
trademarks of Dental Wings.

29-0401-EN v1.0

© 2020 Dental Wings Inc.
A Straumann Group Brand.

T 514 807.8485 | 1 888 856.6997 (toll free in U.S. & Canada)
E info@dentalwings.com
W dentalwings.com/brochures

