

# PC specifications to run DWOS software

Version 2.5      15/04/2020

---

## Recommended PC specifications to run DWOS® Software

- CPU: Quad Core 3.30 GHz [Intel Core i5 3.30 GHz]
- Memory RAM: 16 Gigabytes
- Graphic card: ATI Radeon HD6570, 1GB
- Hard drive (HD): 500 gigabytes 7200 rpm
- Network with internet connection
- Monitor 22" wide 1080p (recommended)
- USB keyboard
- USB mouse (two buttons and a wheel)

## Operating system

- Windows 7 (64-bit edition)
- Windows 10 (64-bit edition)

## Optional components

- 3Dconnexion navigation devices